

Pevco

TEC-6™ Trackable Ergonomic Carrier, 6"




Exclusively for 6" Healthcare
Pneumatic Tube Systems

Pevco TEC-6™ Trackable Ergonomic Carrier, 6”

Exclusively for 6” Healthcare Pneumatic Tube Delivery Systems

Pevco TEC-6's ergonomic design incorporates well-placed, comfortable hand-holds and contoured ends so healthcare staff can hold and load carriers with ease. Red/green latch status indicators help guarantee that carriers are closed properly.

Pevco TEC-6™ carriers are durable, maintain lasting clarity and are leak resistant. With long-wearing nylon glide bands, tough hardened steel hinges and Pevco's patented leak-resistant seal, the TEC-6 is built to stand up to the most demanding 24/7 hospital environment.



Combine Pevco TEC-6 with Pevco Link™ and benefit from a real-time audit trail for every delivery.

Know what was sent, when it was sent, when it arrived and where it arrived – and even who sent and received it with barcode tracking. Transaction data is easily accessed from a station touch control and more detailed information is available through an accompanying web application.



Foam inserts are available to protect fragile items.

TEC-6

Long-wearing nylon glide bands




Dependable low-profile latches

Pevco TEC-6™ also includes integrated barcodes to enable real-time tracking of every delivery.



When combined with Pevco Link™, the contents of every TEC-6 can be tracked from its origin to destination providing a real-time audit trail.



Built to stand up to the most demanding hospital environment with long-lasting, tough hardened steel hinges and Pevco's patented leak resistant seal



Available in three styles, foam inserts help protect fragile materials



Available in three carrier color options shown below:



Conveniently placed hand holds

Red/green latch status indicators

The contoured end allows the carrier to stand without support, and makes them easy to grab.

TEC-6

Trackable

Ergonomic

Carrier, 6"

Pevco™

1401 Tangier Drive, Baltimore, Maryland 21220
Tel: 410.931.8800 800.296.7382 Fax: 410.931.4660

www.pevco.com

