

Pevco TEC-6™

TRACKABLE ERGONOMIC CARRIERS



Pevco TEC-6™ Trackable Ergonomic Carriers

Exclusively designed for 6" pneumatic tube delivery systems.



Durable

Built to stand up to demanding hospital environments, with tough, hardened steel hinges and long-wearing nylon glide bands.



Leak-resistant

Made with a leak-resistant seal to help avoid spills, costly downtime and cleanup.



Ergonomic

Designed with comfortable hand-holds, low profile rotating latches and belled ends that allow carriers to stand on end.



Trackable

Includes integrated RFID tags and barcodes for real-time tracking to improve compliance and streamline empty carrier allocation.
(For Pevco Atlas™ system with Pevco Auto-ID™)



Antistatic

Long-wearing nylon glide bands ensure smooth transport, and include copper thread to eliminate static build up.

Call **800.296.7382**

for more information or to order.

With the addition of Pevco Auto-ID™


TEC-6 RFID tags enable enhanced delivery monitoring and streamlined carrier inventory management.

When combined with Pevco Link™

Contents of every TEC-6 can be tracked from origin to destination, providing an end-to-end audit trail.

For systems with Pevco Pilot™

Managers can create rules that govern how deliveries are managed, reducing turnaround times and errors.



Protect Fragile Items

Foam inserts are available in egg crate, 1" thick and 1/2" thick



Color Options

Available in black, yellow, blue, red and green





pevco.com

1401 Tangier Drive
Baltimore, Maryland 21220

Tel

410.931.8800

800.296.7382

Fax

410.931.4660